

Notification to consumers and HCPs: LUCIRA® by Pfizer COVID-19 & Flu Home Test and Lucira® by Pfizer for COVID-19 & Flu Test (For use under Emergency Use Authorization (EUA) only)

The notification is to inform consumers and HCPs that certain lots of the test kits include nasal swabs that expired before the labeled expiration date printed on the outer test package. The impacted kits expire between February and June 2025, while the nasal swab expired in January 2025.

Impacted Test Kit Lots:

Name	Kit Lot #
LUCIRA by Pfizer for COVID-19 and Flu Test	K10A011512233M4
LUCIRA by Pfizer COVID-19 & Flu Home Test	K10A112012233M1 K10A112012233M2 K10A112012233M3 K10A112012233M4 K10A112012233M5 K10A112012233M6 K10A112012233M7 K10A112012233M8

The affected products, LUCIRA® by Pfizer COVID-19 & Flu Home Test and Lucira® by Pfizer for COVID-19 & Flu Test® are single-use home test kits for detecting SARS-CoV-2, influenza A, and influenza B. These test package includes a detection module, sample preparation tube, and one package with a sterile swab manufactured by Copan Industries, Inc.

- Please check the expiration date on the swab packing (See Figure 1, below).
- If the swab in your kit is expired, please do not use and call LUCIRA® by Pfizer Customer Support [1-888-LUCIRA-4 \(582-4724\)](tel:1-888-LUCIRA-4). Monday-Sunday 8:00 am-8:00 pm ET to receive a replacement kit

Figure 1. The expiry date on the Copan swab package is highlighted in red. The expiry date and lot number for the Lucira by Pfizer kit is highlighted in blue.

